



PSU NMR Facility Acquiring ^1H NMR spectrum (Quick instructions for Topspin)

If you need assistance, please contact:

- Dr. Carlos Pacheco -- Room#14 - cnp14@psu.edu; cpacheco223@gmail.com;
Ph# office: (814)863-1182; cell: (609) 240-5957.
- Dr. Debashish Sahu -- Room#10 - desahu@psu.edu

Place your sample in the spinner and carefully adjust the NMR tube depth in the gauge.

- 1) **Turn on** the lift air using BSMS (or type *ej* <enter> at the command line).
 - a. If BSMS window is missing, type *bsmsdisp* <enter> at the command line.
 - b. ***Always wait for the sound of the lift air at full power before placing your sample on top of the magnet!!!***
- 2) **Turn off** lift air from BSMS (or type *ij* <enter> at the command line).
- 3) Create a new data set by typing *new* <enter> at the command line or click  on the toolbar (or type *edc* at the command line)



- NAME -- (filename for your compound)
- EXPNO -- suggestion: Set 1 for a ^1H , 2 for a ^{13}C
- PROCNO -- set to 1
- DIR -- set to /opt/topspin/
- USER -- should be <your_Professor-name>
- SOLVENT -- your deuterated solvent from the pull-down menu.
- EXPERIMENT -- "Use current params"
- TITLE -- the title of your sample (can be as descriptive as possible)

- 4) Read a good shimfile with **rsh**<enter>: choose the **stdshims** file from the list, which is the standard shimfile name for all NMRs.
- 5) Type **lock** <enter> at the command line and choose the appropriate solvent from the lock table.
 - a. If the lock window is missing, type **lockdisp** <enter> at the command line.
- 6) After the lock process is finished, **shim z1 and z2**.
- 7) Type **rpar PROTON all** <enter> at the command line.
- 8) Type **getprosol** <enter> at the command line.
 - a. For instruments running **Xwinnmr**, type **eda** <enter> and set prosol to **TRUE**.
- 9) (OPTIONAL) Adjust **SW** and **o1p** values.
- 10) Type **rga** <enter> at the command line
 - a. Adjust the number of scans - **NS** - as appropriate if necessary; default for ¹H is 16, for ¹³C is 1024).
- 11) Type **zg** <enter> at the command line to start the experiment.
- 12) Type **acqu** <enter> at the command line if you do not see your FID (raw NMR data).
- 13) Your NMR data is automatically saved at the end of the experiment.

DATA PROCESSING

((It is highly recommended to carry out the NMR processing on the *off-line* PC at the NMR Facility – room #10, where you can also back your NMR data up))

- 14) Type **ef; apk** <enter> at the command line to perform Apodization (**e**), Fourier (**f**) Transformation, and automatic phase adjustment (**apk**).
- 15) If needed, further optimize the phase of the spectrum manually, **using the phase correction icon**.


Click on the phase icon:





Left-click **0**, hold the mouse button to correct phase at the vertical cursor.

Left-click **1**, hold the mouse button to correct phase away from the cursor.


Click on the SAVE icon to exit the phasing window.



16) Calibrate your spectrum using a known chemical shift (residual solvent peak or TMS). Click on the icon  at the toolbar. Set the reference and return.


17) Type **abs** to flatten the baseline. There is a manual baseline flattening: 

18) Integrate the spectrum manually with  -- This takes you into a subroutine:



Integrate the spectrum, save the integrals before exit the window.

19) Peak picking can be done with the icon  and drawing a green box around the peaks:



Save the peaks before exiting the window (using the same icon as the phasing and integration).

20) Type ***plot*** <enter> to open the plot editor (for instruments running TopSpin) and then print.

- a. For instruments running Xwinnmr type ***view*** <enter> to print preview and type ***plot*** <enter>.

21) Print.

22) Save your data on the PC at Room #10

- a. Save your NMR data **ALSO** on a thumb drive.
 - i. Email the NMR data is also possible at the off-line PC.