This handout was created for stroke survivors and their family members and caregiver to provide easy and simple instructions on what equipment they might need during discharge in acute care setting hospital. The following adaptive equipment were recommended based on the experience in the hospital. The therapists can print this information and hand it to stroke survivors with check mars if necessary.

L Tub guard bathtub safety rail

Tub-guard bathtub safety rails offer you to move around safely inside their tub. This will prevent falling during tub transfer. DO NOT PULL on rail, it is for pushing off of.



Hand-held shower with ergo handles

The long hand-held shower provides a comfortable,

nonslip grip. This equipment is helpful during shower, where you plan to take a shower while sitting up on a shower chair or tub transfer chair. This equipment will be easier to reach parts of



your body and able to take shower independently.

Leg lifter Long handled Reacher This equipment will allow to pick up objects without The leg lifter strap allows you to lift up your leg when bending; which will prevent falling. Also, this can help to your leg muscles are not able to do so on their own. dress or undress pants. "Clip on" Reacher is recommended for you if you are using walkers. Contoured long Transfer bench with suction cups If you need extra support getting in and out of the bathtub handle sponges or shower, transfer bench will be helpful for you. The bath Designed to promote a more independent those with limited chair is reversible with height-adjustable legs to fit most reach. Handled Sponge comes with bathtubs. There are suction cups under the two legs inside a sponge head that gently washes tub for and scrubs all the hard to reach stabilizing areas. This equipment will give you the bench. a way to shower independently. Grab Bars Provides a reliable handhold while showering, bathing, or using the toilet.

Non-slid rubber mats with suction cups underneath



Provides effective slip resistance in the bathtub; suction cups for reliable non-slid traction.



Shower chair /stool

Allows you to prevent fall while showering. Adjustable legs allow to customize chair based on specific preference and drainage holes in bench reduce slipping. Also saves energy or reduce stress on lower extremities. Shower chair with back is recommended if you need back support.

Images retrieved from:

https://bestchoicesforseniors.com/best-non-slip-bath-bathroom-mats-rugs-for-elderly-seniors/

https://thehelpinghome.com/aging-in-place-at-home-definitive-guide-part6/

http://getgrabbarsinstalled.com/reno

https://www.completecareshop.co.uk/bathing-aids/sponges-and-washing-aids/

https://www.arthritissupplies.com/leg-lifter.html

 $https://www.allegromedical.com/daily-living-aids-c519/aluminum-reacher-with-magnetic-tip-p191918.html?CS_003=9164468\&CS_010=8ab281020bb66dff010bb67e52885804\&gclid=EAIaIQobChMI58bymcub4AIVEdVkCh2nDQpuEAQYASABEgIBqvD_BwE#DMI-1764$

 $https://express.google.com/u/0/product/14329802860269338840_0_5826996?utm_source=google_shopping\&utm_medium=tu_cu\&utm_content=eid-product/source=google_shopping\&utm_medium=tu_cu\&utm_content=eid-product/source=google_shopping\&utm_medium=tu_cu\&utm_content=eid-product/source=google_shopping\&utm_medium=tu_cu\&utm_content=eid-product/source=google_shopping\&utm_source=google$

lsjeuxoeqt>im=CP6Bw9L5rpPmQhDtmILC-

reA328YwOrWFiIDVVNEKKC6mOMFMLTT4wI&utm_campaign=5826996&gclid=EAIaIQobChMIu9ibosub4AIVC8RkCh3oSw45EAQYASABEgKbrPD BwE

This handout was created in partial fulfilment of OTD degree on February 1, 2019.



This handout was created in partial fulfilment of OTD degree on February 1, 2019.